



FIGURE S1. IMMUNODETECTION IN CASHEW



- 1- Untreated
- 2- Boiled 30 min
- 3- Boiled 60 min
- 4- Auto 121°C 15 min
- 5- Auto 121°C 30 min
- 6- Auto 138°C 15 min
- 7- Auto 138°C 30 min

FIGURE S2. IMMUNODETECTION IN PISTACHIO




Figure S3. Western blot de hazelnut (3 patients)





Figure S4. Western blot de almond (pool 4 patients)

Autoclave



Altas presiones





FIGURE S5. IMMUNODETECTION IN CHESNUT



IgE ELISA