

Thermo fisher

BenchMark Pre-stained Protein Ladder
Cat. no. 10748-010




Figure 4C SESN 1 (57kd) lane 2, 3 and 4

180kDa
115 kDa
82 kDa
64 kDa
49 kDa
37 kDa
26 kDa
19 kDa
15 kDa
6 kDa




Figure 4C SESN 2 lane 2, 3 and 4

180kDa
115 kDa
82 kDa
64 kDa
49 kDa
37 kDa
26 kDa
19 kDa
15 kDa
6 kDa

Figure 4C B actin lane 2, 3 and 4

NeoPRO 10 Prestained Protein Ladder NB-59-0002





Figure 5B SESN1 lane 2 and 3




Figure 5B SESN2 lane 2 and 3




Figure 5B actin lane 1 and 2




Figure 7A P53 lane 2 and 3




Figure 7A p21 lane 2 and 3




Figure 7A B actin lane 2 and 3




Figure 7B P53 lane 2 and 3




Figure 7B p21 lane 2 and 3




Figure 7C B actin lane 4 and 5




Figure 7C MMP13 lane 2 and 3




Figure 7C B actin lane 2 and 3