

Article

RgIA4 Prevention of Acute Oxaliplatin-Induced Cold Allodynia Requires $\alpha 9$ -Containing Nicotinic Acetylcholine Receptors and CD3⁺ T-Cells

Peter N. Huynh ^{1,*}, Sean B. Christensen ¹ and J. Michael McIntosh ^{1,2,3}

¹ School of Biological Sciences, University of Utah, Salt Lake City, UT 84112, USA

² George E. Wahlen Veterans Affairs Medical Center, Salt Lake City, UT 84112, USA

³ Department of Psychiatry, University of Utah, Salt Lake City, UT 84112, USA

* Correspondence: peter.huynh@utah.edu

Supplemental Information

Table S1: Commercial antibodies used in flow cytometry analysis and CD3⁺ cell depletion accompanied with test concentrations.

Flow Cytometry Antibodies								
Target	Fluorophore	Clone	$\mu\text{g}/\text{test}$	Host	Reactivity	Isotype	Manufacturer	Catalog No.
Viability	FVS575V	N/A	0.59	N/A	N/A	N/A	BD Horizon	565694
CD45	AlexaFluor 594	30-F11	0.25	Rat	Mouse	IgG2c	Bio-Rad	128011
CD11b	PE/Cyanine7	M1/70	0.625	Rat	Mouse	IgG2b	Invitrogen	25-01122-81
CD3	APC	KT3	10	Rat	Mouse	IgG2a	Bio-Rad	MCA500APC
CD4	Pacific Blue	RM4-5	10	Rat	Mouse	IgG2a	Bio-Rad	MCA2691PB
CD8	PE-Cy5.5	53.6.7	0.625	Rat	Mouse	IgG2a κ	ThermoFisher	35-0081-80
CD25	PE	PC61.5.3	5	Rat	Mouse	IgG1	ThermoFisher	MA5-17817
In vivo depletion antibodies								
Target	Fluorophore	Clone	mg/kg	Host	Reactivity	Isotype	Manufacturer	Catalog No.
CD3 ϵ	-	145-2C11	1000	Armenian Hamster	Mouse	IgG1	BioXCell	BE0001-1