

Supplementary Materials: The Role of Sarcosine, Uracil, and Kynurenic Acid Metabolism in Urine for Diagnosis and Progression Monitoring of Prostate Cancer

Georgios Gkotsos, Christina Virgiliou, Ioanna Lagoudaki, Chrysanthi Sardeli, Nikolaos Raikos, Georgios Theodoridis and Georgios Dimitriadis

Data regarding AUCs, SEs, p and CIs from Receiver operating characteristic analysis are presented in Tables S1–6

Table S1				
Endogenous metabolites	AUC	SE	p	95% CI
Sarcosine	0,47	0,052	0,554	0,37-0,57
Kynurenic acid	0,44	0,05	0,251	0,34-0,54
Uracil	0,59	0,05	0,066	0,49-0,69

Table S2				
Endogenous metabolites	AUC	SE	p	95% CI
Sarcosine	0,38	0,057	0,044	0,27-0,49
Kynurenic acid	0,41	0,057	0,128	0,31-0,52
Uracil	0,47	0,059	0,594	0,35-0,58

Table S3				
Endogenous metabolites	AUC	SE	p	95% CI
Sarcosine	0,46	0,059	0,473	0,34-0,57
Kynurenic acid	0,62	0,062	0,041	0,49-0,74
Uracil	0,54	0,059	0,492	0,42-0,66


Table S4				
Endogenous metabolites	AUC	SE	p	95% CI
Sarcosine	0,48	0,064	0,819	0,36-0,61
Kynurenic acid	0,51	0,064	0,858	0,38-0,63
Uracil	0,54	0,064	0,525	0,42-0,67

Table S5				
Endogenous metabolites	AUC	SE	<i>p</i>	95% CI
Sarcosine	0,51	0,081	0,927	0,35-0,67
Kynurenic acid	0,5	0,083	0,985	0,34-0,66
Uracil	0,54	0,082	0,595	0,38-0,71


Table S6				
Endogenous metabolites	AUC	SE	<i>p</i>	95% CI
Sarcosine	0,52	0,081	0,84	0,36-0,68
Kynurenic acid	0,57	0,08	0,36	0,42-0,73
Uracil	0,52	0,083	0,777	0,36-0,69

Supplementary Figure S1

a **Sarcosine**



b **Uracil**



c **Kynurenic acid**




Figure S1: Chromatographic peaks of a) sarcosine, b)uracil and c)kynurenic acid in urine samples.