

## Supporting Information

# Analytical Method and Stability Study for Oral Suspension of Idebenone in Syrspend

Emanuele Porru <sup>1</sup>, Filomena Piro <sup>2</sup>, Rossana Comito <sup>1</sup>, Anastasiia Mosendz <sup>3</sup>, Elirosa Minniti <sup>2</sup>, Matteo Conti <sup>4</sup>, Alessandra Stancari <sup>2</sup> and Francesco Saverio Violante <sup>1,\*</sup>

<sup>1</sup> Occupational Medicine Unit, Department of Medical and Surgical Sciences, Alma Mater Studiorum—University of Bologna, 40138 Bologna, Italy; emanuele.porru2@unibo.it (E.P.); rossana.comito2@unibo.it (R.C.)

<sup>2</sup> Clinical Pharmacy Production and Research, IRCCS Azienda Ospedaliero—Universitaria di Bologna, 40138 Bologna, Italy; filomena.piro@aosp.bo.it (F.P.); elirosa.minniti@aosp.bo.it (E.M.); alessandra.stancari@aosp.bo.it (A.S.)

<sup>3</sup> Department of Analytical Chemistry, Taras Shevchenko National University of Kyiv, 01033 Kiev, Ukraine; anastasiia.mosendz@studio.unibo.it


<sup>4</sup> Local Unit of Imola, Department of Public Health, Health Service of the Emilia—Romagna Region, 40026 Imola, Italy; matteo.conti@ausl.imola.bo.it

\* Correspondence: francesco.violante@unibo.it; Tel.: +39-051-2142611; Fax: +39-051-2142609

**Table S1. Idebenone mean concentration value expressed as relative percent of recovery after forced degradation.**

Stress condition	Duration	IDB concentration (%)
Acid hydrolysis	10 days	100%
Base hydrolysis	24 hours	75.2%
Oxidation	7 days	90.3%
Thermal (50°C)	10 days	100 %
Photolytic	11 days	100 %

**Figure S1. Calibration curve in Mobile phase (A) and in matrix (B)**



**Table S2. Accuracy (bias %) and precision (CV%) intra- (n=3) and inter-day (n=9) repeatability at three concentrations (low, medium, high).**

Compound	Intra-day (n=3)						Inter-day (n=9)					
	QClow		QCmedium		QChigh		QClow		QCmedium		QChigh	
	Bias	CV	Bias	CV	Bias	CV	Bias	CV	Bias	CV	Bias	CV
	%	%	%	%	%	%	%	%	%	%	%	%
IDB	5,2	4,3	4,5	2,3	4,5	1,7	5,2	2,8	4,2	3,6	4,8	3,8